


## TT30-D: Tub Temperature Testing Card, Stock Dolphin

TT30-D: front



TT30-D: back

- Make sure hot water heater is set at 120 degrees.
- Fill tub with hot and cold water at the same time.
- Fill tub no higher than baby's belly button.
- Use this card and your hand when checking water temperature. If skin flushes, water is too hot.
- Put your child in the water with their back to the faucet so they cannot turn the water on.
- Never leave baby unattended.
- Empty bath as soon as you are done.

**This card is a guide and is to be used in conjunction with the hand method**

Please remove all guideline indicators before submitting artwork



= Final Cut Size: 2.125" x 3.375"



= Ad Copy Area: 1.875" x .375"  
Art must be within clearances - 0.125" on all sides.